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


Plasma Drug Concentration Time Curve

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Anti-idiotypic Antibodies

Supporting Immunogenicity and Pharmacokinetics Analysis

Method Validated and Protocol Offered

- High Affinity
- High Specificity
- High Stability

カタログ番号	抗原	中和活性	親和性 (Ko.nM)	用途
TRB-Y1b	Trastuz*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA with TRB-Y5b Neutralizing assay Indirect ELISA
TRB-Y5b	Trastuz*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA with TRB-Y1
ADB-Y19	Adalim*mab F(ab') ₂	Neutralizing Antibody	0.0013	ADA assay Neutralizing assay Indirect ELISA
ADB-Y23b	Adalim*mab F(ab') ₂	Non-Neutralizing Antibody	/	PK bridging ELISA with ADB-BY17 Indirect ELISA
ADB-BY17	Adalim*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA with ADB-Y23
CEB-Y27	Cetux*mab F(ab') ₂	Neutralizing Antibody	0.007	ADA assay Neutralizing assay Indirect ELISA
CEB-Y28	Cetux*mab F(ab') ₂	Neutralizing Antibody	0.0015	ADA assay Neutralizing assay Indirect ELISA
CEB-Y29	Cetux*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA with CEB-BY31 Neutralizing assay Indirect ELISA
CEB-Y31	Cetux*mab F(ab') ₂	Non-Neutralizing Antibody	0.421	ADA assay Indirect ELISA
CEB-BY31	Cetux*mab F(ab') ₂	Non-Neutralizing Antibody	/	PK bridging ELISA with CEB-Y29
RIB-Y36	Ritux*mab F(ab') ₂	Neutralizing Antibody	0.01	ADA assay Neutralizing assay Indirect ELISA
RIB-Y37	Ritux*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA Neutralizing assay Indirect ELISA
BEB-Y12	Bevaciz*mab F(ab') ₂	Neutralizing Antibody	0.0828	Neutralizing assay Indirect ELISA
BEB-Y9	Bevaciz*mab F(ab') ₂	Neutralizing Antibody	1.92	ADA assay Indirect ELISA
BEB-Y10	Bevaciz*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA with BEB-BY13 Neutralizing assay Indirect ELISA
BEB-Y13	Bevaciz*mab F(ab') ₂	Neutralizing Antibody	/	PK bridging ELISA with BEB-Y10

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カタログ番号	製品説明	標準
EPH-V1	ELISA Assay Kit for Anti-PD-1 h-mAb in Human Serum	96/480 tests
EPM-V1	ELISA Assay Kit for Anti-PD-1 h-mAb in Mouse Serum	96/480 tests
EPC-V1	ELISA Assay Kit for Anti-PD-1 h-mAb in Monkey Serum	96/480 tests
EHH-V1	ELISA Assay Kit for Anti-HER-2 h-mAb in Human Serum	480 tests
EHM-V1	ELISA Assay Kit for Anti-HER-2 h-mAb in Mouse Serum	480 tests
EHC-V1	ELISA Assay Kit for Anti-HER-2 h-mAb in Monkey Serum	480 tests
ECH-V1	ELISA Assay Kit for Anti-CTLA-4 h-mAb in Human Serum	96/480 tests
ECM-V1	ELISA Assay Kit for Anti-CTLA-4 h-mAb in Mouse Serum	96/480 tests
ECC-V1	ELISA Assay Kit for Anti-CTLA-4 h-mAb in Monkey Serum	96/480 tests

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






用途	発現方法	概略図	長所	制限事項
PK assay	Development of an anti-idiotypic antibody targeting the variable region of the drug	Secondary antibody Drug Anti-idiotypic antibody	Only one anti-idiotypic antibody needs to be developed; Low costs	Specificity is general, only suitable for drugs containing Fc domains
	Development of a pair of anti-idiotypic antibodies targeting the variable region of a drug	Anti-idiotypic antibody Anti-idiotypic antibody Drug	High specificity and accuracy	Development is difficult
ADA assay	Development of a multi-antibody targeting the drug's variable region	HPR-labeled drug Anti-idiotypic antibody Drug	/	/

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用途	標的	発現方法	概略図
PK assay	Complete bispecific antibody	Development of paired anti-idiotypic antibodies for two units	 <p>Bispecific antibody Anti-idiotypic antibody HRP-labeled anti-idiotypic antibody</p>
	A unit	Development of paired anti-idiotypic antibodies, targeting the A unit	 <p>HRP-labeled anti-idiotypic antibody Bispecific antibody Anti-idiotypic antibody</p>
	B unit	Development of paired anti-idiotypic antibodies, targeting the B unit	 <p>HRP-labeled anti-idiotypic antibody Bispecific antibody Anti-idiotypic antibody</p>
ADA assay	Complete bispecific antibody	Development of a multi-antibody against full-length antibody (variable region)	 <p>HRP-labeled drugs Anti-idiotypic antibody Bispecific antibody</p>

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用途	標的	発現方法	概略図
PK assay	Antibody-drug conjugate	Development of paired anti-idiotypic antibodies for antibodies and small molecules	
	Antibody	Development of paired anti-idiotypic antibodies that target antibodies	
	Free small molecules	/	LC-MS
ADA assay	Antibody-drug conjugate	Develop a multi-antibody for ADC	

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用途	標的	発現方法	方法
PK Assay	scFv	Development of an anti-idiotypic antibody for scFv	Cell-Based Assay
ADA Assay	scFv	Developed a multi-antibody for scFv	Cell-Based Assay/ELISA

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	抗原の補足の種類	抗イディオタイプ抗体の補足の種類	サンドイッチ酵素結合免疫吸着測定法
Schematic diagram of detection	<p>HRP/Biotin-Labeled Anti-Fc secondary antibody</p> <p>Antibody drug</p> <p>Target protein</p>	<p>Anti-Fc secondary antibody</p> <p>Antibody drug</p> <p>Anti-idiotypic antibody</p>	<p>Anti-idiotypic antibody</p> <p>Antibody drug</p> <p>Anti-idiotypic antibody</p>
Solid phase coating	Target protein	Anti-idiotypic antibody	Anti-idiotypic antibody
Target	Antibody drug	Antibody drug	Antibody drug
Antibody detection	Anti-Fc secondary antibody	Anti-Fc secondary antibody	Anti-idiotypic antibody
Advantages	No need for additional development of anti-idiotypic antibodies. Simple and time-saving protocol	Only one anti-idiotypic antibody needs to be developed. Saving costs	High accuracy. Suitable for all types of biological drugs
Limitations	Low degree of stability and certainty	Average specificity, susceptible to matrix effects, only suitable for antibody drugs with Fc region	High development difficulty

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Testing method	Coated	Sample	Testing
Antigen capture ELISA	CD20	Antibodies to be tested	Goat anti-human IgG
Anti-idiotypic capture ELISA	Anti-Ritux [®] mab Antibodies	Antibodies to be tested	Goat anti-human IgG
Bridging ELISA by anti-idiotypic antibodies	Anti-Ritux [®] mab Antibodies	Antibodies to be tested	Biotinylated Anti-Ritux [®] mab Antibodies


Figure1. Detection of Ritux[®]mab by antigen-capture ELISA (0.1% serum).


Figure2. Detection of Ritux[®]mab by anti-idiotypic capture ELISA (0.1% serum).


Figure3. Detection of Ritux[®]mab by anti-idiotypic bridging ELISA (10% serum).

Testing method	Linear range (µg/mL)	Sensitivity (µg/mL)	Advantage	Disadvantage
Antigen-capture ELISA	—	—	Simple method and good versatility	High background, no activity
Anti-idiotypic capture ELISA	0.156-10	0.156	Solve the difficulty in obtaining CD20, simple method	High background, only suitable for Rituxa biosimilar
Bridging ELISA by anti-idiotypic antibodies	0.012-0.78	0.012	Solve the difficulty in obtaining CD20, good sensitivity and low background	Only applicable to Rituxan biosimilar

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	抗体医薬品の補足の種類	ブリッジングの種類
Schematic diagram of detection		
Solid phase coating	Antibody drug	Antibody drug
Target	Anti-drug antibody	Anti-drug antibody
Antibody detection	Anti-Fc secondary antibody	Antibody drug
Advantage	Easy to develop	High accuracy
Limitations	Need to replace the detection antibody, accuracy is generally affected by the matrix. For testing animal serum only	Narrow linear range. Difficult to develop

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Anti-Adalimumab Antibodies
(ADB-Y19)


Figure 1. Anti-Adalimumab Antibodies bridging ELISA for Anti-Drug Antibody (ADA) assay development. Immobilized adalimumab at 1 µg/ml, added increasing concentrations of Anti-Adalimumab Antibodies (Cat. No. ADB-Y19, 10% human serum) and then added biotinylated adalimumab at 5 µg/ml. Detection was performed using HRP-conjugated streptavidin with a sensitivity of 0.6 ng/mL.

Anti-Rituximab Antibodies
(RIB-Y35)


Figure 2. Anti-Rituximab Antibodies bridging ELISA for Anti-Drug Antibody (ADA) assay development. Immobilized rituximab at 5 µg/ml, added increasing concentrations of Anti-Rituximab Antibodies (Cat. No. RIB-Y35, 10% human serum) and then added biotinylated rituximab at 5 µg/ml. Detection was performed using HRP-conjugated streptavidin with a sensitivity of 20 ng/mL.

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AS-CN03E01-OPE-20210702

*Where Proteins and Innovation
Advance Biomedicine*

Copyright Statement

“

This material is copyrighted by the Company. All rights in this material are reserved by the Company. Unless otherwise indicated in writing, all material in this material is copyrighted by the Company. No part of this material may be copied, photocopied or reproduced in any form or redistributed to any other person or used in any other manner which infringes the Company's copyright without the prior written authorisation of the Company.

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